

# HELPING YOUR DOG STAY CALM

## An anxious dog can't help the way they feel

In dogs with separation anxiety, the area of the brain responsible for processing emotions and detecting fear goes into overdrive when they can't be with their owners. It's the equivalent of a full-blown panic attack in a human.

## Calmer Canine restores emotional balance to a brain in turmoil


Calmer Canine's treatment is invisible, sensation-free, and very safe. Just two 15-minute treatments are needed each day for 4 to 6 weeks. Your dog simply wears their vest and device and goes about their day. Treatment is given while you're home, not while you're away.

# WHERE TECHNOLOGY MEETS LOVE



Science has proven that a powerful link exists between anxiety and inflammation in the brain. An anxious brain is out of balance—with overactive brain cells—producing harmful substances causing inflammation and cell death.

Utilizing a specifically tuned microcurrent signal delivered by the device placed over the dog's head, Calmer Canine stimulates the brain's own ability to reduce inflammation and produce beneficial "feel good" chemicals (like endorphins)—restoring natural balance, bringing your dog back to a calm state.

Calmer Canine's treatment, based on FDA-cleared technology currently in use to treat people with anxiety, depression, and ADHD, is much more desirable than other treatments that don't have lasting effects, including medications with unwanted "dopey dog" side effects.



Calmer Canine helps anxious dogs regain and maintain calm with just two 15-minute at-home treatments a day.



Calmer Canine helps reduce inflammation, boost production of "feel good" chemicals, and restore natural balance to a dog's brain.

# MEET THE CALMER CANINE SYSTEM



- 1 BATTERY-POWERED DEVICE**  
The device features an applicator halo, which emits the microcurrent signals. The applicator halo is available in two diameters (5 inches and 7 inches).
- 2 EASY-TO-USE VEST**  
The comfortable, washable vest holds the device in the ideal position during treatment—allowing for hands-free operation. The device is available with or without the vest.
- 3 SIZES FOR EVERY DOG**  
To ensure the best fit, the vest comes in five sizes: XS, S, M, L, and XL.

- 4 OR USE WITHOUT THE VEST**  
If your dog dislikes wearing garments or they're bothered by the sound of hook-and-loop fasteners, you can hold the device in place during treatments.

